

Tip infection control procedure

After every patient

It is of utmost importance to clean the scanner tip after every patient to ensure proper infection control.

1. **Wipe scanner body down WITH TIP ATTACHED to avoid any germicide coming into contact with the sensitive connection points or lenses.**
2. Remove the scanner tip from the scanner. Place some protective film or similar cover over the scanner body's connection points and lens to prevent ingress of liquids or dust when a tip isn't on the scanner.
3. Soak the tip in warm tap water (40° - 50°C [104° - 122°F]) for 5 minutes.
4. Using a soft drinking straw/bottle brush, rinse and brush the **external surface** of the tip under warm running water (3-3.5L/min @ 30°-35°C [86°-95°F]) for a minimum of 15 seconds, with emphasis on the tip alignment indicator.

A typical water faucet should be able to produce a flow rate of between 4-8 Lpm.

5. Rinse and brush around the **edge and surrounding crevices of the mirror** under running water for a minimum of 30 seconds.
To avoid scratching the mirror, be careful that only the bristles touch its surface.



The example picture is shown without running water for visibility purposes. Use the bristles of the small brush on all crevice edges of the mirror.

6. After brushing, rinse the mirror and surrounding crevice under running water for additional 10 seconds.
7. Wipe the mirror with an optical lens wipe.
 - If there is any visible debris remaining after performing the cleaning steps above repeat the process.
 - If there is any debris that cannot be removed or visible deterioration such as cracks or discoloration, the tip should be disposed of and replaced. Contact your local distributor to purchase replacement.
8. Dry the tip contact area with an optical lens wipe.
9. Place each scanner tip individually into an autoclave pouch.
10. Place one to three pouches per tray or cassette. **DO NOT stack tips on or around other metal instruments.**
11. Select the autoclave cycle with a minimum sterilization temperature of 132° C (269° F). Minimum time and temperature values listed reflect the specific discrete values used to validate the sterilization process for the device.
 - **Use distilled water in the autoclave or tip life may be dramatically reduced as a result.** The high temperatures of the autoclave could cause the chlorine in tap water to corrode the tip contacts.
 - For **Pre-Vac** autoclaves minimum sterilization time is **3 minutes**.
 - For **Gravity** autoclaves minimum sterilization time is **10 minutes**.
 - **Autoclave sterilization cycle times exceeding 10 minutes could be considered excessive and may damage the tip.**
 - For any autoclave, use a minimum dry time of 30 minutes.
12. After cycle completion, store for later use.

Manufacturer's specifications, additional information, and cautions

- Prior to scanning, visually confirm the scanning mirror and scanner connection surfaces are clean and dry.
- To clean either surface, use an optical cloth with a small amount of alcohol to gently remove residue.
- Ensure each surface is dry before connecting it to the scanner.



Autoclave cycle times, parameters, and names differs among manufacturers and autoclave types. Refer to your manufacturer's manual for sterilizing wrapped instruments. **Autoclave sterilization cycle times exceeding 10 minutes could be considered excessive and may damage the tip.** To verify that minimum cycle parameters are met, see the User's manual of the autoclave manufacturer. For more information, refer to the User's manual or call customer support.



DO NOT use a chemical sterilant.
DO NOT place in a ultrasonic cleaner.
DO NOT pouch with other instruments.
DO NOT stack tips on or around other metal instruments.

Excessive cycle times may prematurely reduce the life and function of the tip.

Scan QR code for your language version

English (en)	Deutsch (de)
<p>Tip infection control procedure</p> 	<p>Infektionsschutzverfahren für die Scannerspitze</p> 
Français (fr)	Italiano (it)
<p>Procédure de contrôle des infections de l'embout</p> 	<p>Procedura di controllo infezioni del puntale</p> 
Español (es)	Português (pt)
<p>Procedimiento de control de infecciones de la punta</p> 	<p>Procedimento de controle de infecção de ponta</p> 
Русский (ru)	Polski (pl)
<p>Противомикробная обработка насадок</p> 	<p>Procedura kontroli zakażeń dotycząca końcówki skanera wewnętrznego</p> 
中文 (zh)	日本語 (ja)
<p>取相头感控流程</p> 	<p>チップの感染管理手順</p> 
Suomi (fi)	
<p>Kärjen infektiontorjuntaohje</p> 	

For detailed instructions, please see the user's manual. Rights for changes reserved.

10038775_5/03.2025/en